



Shanghai Yunqing Industrial Co., Ltd.

No.508, Fengyong Road, Fengxian District, Shanghai, China  
Email: yichita@yqgroup.vip

**CAUTION:** The straps on the respirator are completely latex-free. Individuals who are highly sensitive to natural rubber latex can feel free to use.

#### DONNING INSTRUCTIONS:

1. Put the respirator in hand, allowing the head straps to hang freely.
2. Place the respirators against your face, pull the bottom strap over your head, position it around your neck, below the ear.
3. Pull the top strap over your head resting it high at the top of your head.
4. Place the fingers of both hands at the top of nosepiece. Mold the nose area to the shape of the face while moving you finger tips down both sides of the nosepiece.
5. The seal of the respirator on the face should be checked by the wearer prior to entering the work area:
  - a) Cover the front of the respirator with both hands, being careful not to disturb the position of the respirator.
  - b) Inhale sharply, if air leaks around the nose bridge, readjust the nosepiece as described in step 4. if air leaks at the respirator edges, work the straps back along the side of your hand. If you can not achieve proper fit, repeat steps 1-4.
  - c) If no leakage is detected, then work may proceed.

#### FIT-CHECK THE FACE-SEAL AS FOLLOWS:

- Place both hands over the respirator – do not disturb the position of the respirator.
- Exhale for the respirator without exhalation valve or inhale for the respirator with exhalation valve
- If you detect air leak, re-adjust the head straps and/or the nose piece.

The respirator must be used in accordance with OSHA regulations 29 CFR 1910.134 which includes the requirement for fit testing before use.

#### USE INSTRUCTIONS:

1. Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness and may result in sickness or death.
2. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all local government requirements. In the U.S., employers must comply with OSHA 29 CFR 1910.134 which covers medical evaluation, training, and fit testing.
3. Particles that can be dangerous to your health include those not visible to your eyes.
4. Leave the contaminated area immediately and contact supervisor if dizziness, irritation, or other distress occur.
5. Store the respirator away from contaminated areas when not in use.
6. Dispose of used product in accordance with applicable regulations.
7. Inspect the respirator prior to each use to assure that it is in good operation condition. Examine all the respirator parts for signs of damage including the head strap with the welding spots, noseclip, and nose pad. The respirator should be disposed of immediately upon observation of damaged or missing parts. Filtering facepieces are to be inspected prior to each use to assure there are no holes in the breathing zone and no damage has occurred. Holes resulting from ripped or torn filter material around strap welding spots are considered damage. Immediately replace respirator if damaged.

#### USE LIMITATIONS:

1. Do **NOT** use the respirator and enter or stay in a contaminated area under the following circumstances:
  - atmosphere contains less than 19.5% oxygen.
  - atmosphere contains oil aerosols.
  - for protection against gases or vapors.
  - contaminants or their concentrations are unknown or immediately dangerous to life or health.
  - for sandblasting, paint-spray operations, asbestos.
2. Do **NOT** modify or misuse the respirator.
3. Do **NOT** use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good face-seal.

**WARNING:** If you cannot achieve satisfactory fit, do NOT enter the contaminated area. See your Supervisor or contact us via email [at: yichita@yqgroup.vip](mailto:yichita@yqgroup.vip).



Shanghai Yunqing Industrial Co., Ltd.

No.508, Fengyong Road, Fengxian District, Shanghai, China 201414

Email: yichita@yqgroup.vip

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC-	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
		YQD850	
84A-9397	N95	X	ABCJMNOP

#### 1. Protection

N95 - Particulate Filter (95% filter efficiency level)  
Effective against particulate aerosols free of oil,  
time use restrictions may apply

#### 2. Cautions and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.  
B - Not for use in atmospheres immediately dangerous to life or health.  
C - Do not exceed maximum use concentrations established by regulatory standards.  
J - Failure to properly use and maintain this product could result in injury or death.  
M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.  
N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.  
O - Refer to User Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.  
P - NIOSH does not evaluate respirators for use as surgical masks.

**WARNING:** This respirator helps protect against certain contaminants. It may not eliminate the risk of contracting disease or infection. **Misuse may result in sickness or death.**

Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals.

FOR MORE INFORMATION or assistance, please contact  
Shanghai Yunqing Industrial Co., Ltd.